




USER GROUP

To provide comfortable seating to users with mild to moderate postural requirements, to help reduce the formation of pressure ulcers and help to stand or sit users with limited mobility.

APPROPRIATE FOR:



NOT APPROPRIATE FOR:

- Demanding postural requirements
- Users with fixtures and contractures
- Users with a limited extension at the knees
- Users who need to be transported from room to room
- Bariatric users

INTENDED ENVIRONMENT



Domestic use



Care homes or nursing homes



'Reduce Shear' recline to maximise pelvic stability
Integrated self-supply overlay for easy integration of existing cushions*
Compatible with transfer aids including Cricket, ReTurn 7500 and Molift Raiser
Knock down design for easier storage and single person delivery
Pressure reducing fabric on backrest and armrest as standard
Colour matched vapour permeable fabric
Castors for easy manoeuvrability
More depth adjustment options using backrest brackets. Reduces the need for depth adjustment kits
60kg legrest capacity

**For further information on the integrated self-supply overlay please see page 11.*

Name:

Department:

Address:

.....

.....

Telephone:

Email:

Name:

Address:

.....

.....

Telephone:

Email:

User weight:

Pressure Care:

Posture:

Mobility:

Diagnosis:

Prognosis:

Transfer method:

Current seating:





Problems with current seating:

Notes:

.....

.....

.....



SMALL	CR5404 - Duratek	406 mm (16") #1	388 mm (15") #4	457 mm (18")	130 kg
	CR5422 - Beige Vinyl	457 mm (18")	406 mm (16")		
	CR4443 - Black Vinyl	508 mm (20")	431 mm (17") #2		
MEDIUM	CR5406 - Duratek	406 mm (16") #1	439 mm (17") #4	508 mm (20")	160 kg
	CR5424 - Beige Vinyl	457 mm (18")	457 mm (18")		
	CR5444 - Black Vinyl	508 mm (20")	482 mm (19") #2		
LARGE	CR5408 - Duratek	406 mm (16") #1	439 mm (17") #4	559 mm (22")	160 kg
	CR5426 - Beige Vinyl	457 mm (18")	457 mm (18")		
	CR5445 - Black Vinyl	508 mm (20")	482 mm (19") #2		

#1 Depth Adjustment

Kit 2" (50 mm)
CA2416



#2 Height Adjustment

Kit 1" (25 mm)
CA2421



#3 Height Adjustment


Kit 2" (50 mm)
CA2422







#4 Height Reduction

Kit 1" (18 mm)
CA2423






Visco (std)	<p>Prevention: Heat sensitive memory foam moulds to the user providing maximum support and pressure reduction.</p> <p>Suitable for users who are at risk of pressure damage. 160kg (25 stone) weight capacity.</p>	N/A – standard part		
Integrated Self Supply	<p>Allows integration of self-supply cushion of approximate dimensions 17.5" x 17.5" x 3".</p> <p>(Professional input required. See page 11 for instructions).</p>	N/A – standard part		
Cushionair - Alternating Seat Cushion with Pump	<p>powered air cushion for users who have developed pressure damage.</p> <p>Weight capacity: 160 kg.</p>	Fits all sizes CA2400		




The Cushionair cushion is a dynamic pressure relieving cushion which is normally recommended for users who have developed pressure damage*.


The principle is that the cells of the cushion inflate and deflate alternately. This provides pressure redistribution via cyclic changes in loading and unloading













**This is based around anecdotal evidence and Accora's experience; we do not yet have clinical studies to evidence this although these are in progress.*



Minimum (mmHg)	0.00
Maximum (mmHg)	84.57
Average (mmHg)	21.51
Variance (mmHg)	347.87
Standard deviation (mmHg)	18.65
Coefficient of variation (%)	86.99
Standard error (mm)	21.48
Standard error (mm)	28.02
Standard error (mm)	111.74
Regional distribution (%)	100.00





<p>Multi-adjustable Pillow Backrest (std)</p>	<p>Accommodation: accommodates posture that has already developed i.e. kyphosis; our most popular backrest.</p>	<p>N/A Standard Part</p>		
<p>Lateral Supports</p>	<p>used under the Multi-adjustable Pillow Backrest cushions to provide moderate postural support SOLD IN A PAIR</p>	<p>CA2413</p>		<p>QTY 2</p>
<p>Lateral Support Backrest</p>	<p>encourages a good midline sitting position for users whose posture is flexible.</p>	<p>Beige Small - CA2406 Beige Medium - CA2407 Beige Large - CA2408</p>		<p>S M L</p>
<p>Lateral Support Backrest</p>	<p>encourages a good midline sitting position for users whose posture is flexible.</p>	<p>Small - CA2540 Medium - CA2541 Large - CA2542</p>		<p>S M L</p>
<p>Adjustable Lateral Support Backrest</p>	<p>contoured foam backrest with adjustable Velcro Lateral Supports. Provides more aggressive lateral support and mild positioning/correction. Supplied with breathable, pressure reducing fabric.</p>	<p>Small - CA2410 Medium - CA2411 Large - CA2412</p>		<p>S M L</p>
<p>Profiled Headrest</p>	<p>Suitable for users with reduced lateral head control. Fits over the shoulders.</p>	<p>Fits all sizes Beige - CA2419</p>		
<p>Profiled Headrest</p>	<p>Suitable for users with reduced lateral head control. Fits over the shoulders.</p>	<p>Black - CA2530</p>		
<p>2-Button Handset</p>	<p>Suitable for Cognitive Impaired person e.g. dementia (Chair comes Standard w/ 4-button Handset).</p>	<p>Fits all Chair Sizes SC2316</p>		
<p>Legrest Channel</p>	<p>Provides side support on legrest to reduce abduction</p>	<p>Small - CA2520 Medium - CA2521 Large - CA2522</p>		<p>S M L</p>

continued overpage

Right-hand Side
Configura Handset
Holder

A great addition for users with dementia, or similar conditions, who may have difficulties locating the handset. Left or right-hand side option.

CA2556



Left-hand Side
Configura Handset
Holder

Ideal for users with limited movement on one side of their body from conditions like brain injury or stroke. Users can easily reach the handset from any position, including recline and tilt-in-space, while maintaining optimal posture.

CA2554



Lefthand side
Dropdown Armrest

Enables user to do a slideboard transfer in and out of the chair with ease. Can be bolted on at any time if the users needs change.

Duratek - CA2484
Beige Vinyl - CA2492
Black Vinyl - CA2534

Duratek - CA2488
Beige Vinyl - CA2496
Black Vinyl- CA2535



S
M
L

Righthand side
Dropdown Armrest











Enables user to do a slideboard transfer in and out of the chair with ease. Can be bolted on at any time if the users needs change.

Duratek - CA2486
Beige Vinyl - CA2494
Black Vinyl - CA2536

Duratek - CA2490
Beige Vinyl - CA2498
Black Viny - CA2537



S
M
L

Armrest overhang on outside (std)	Provides a good grip to assist the user rolling from side to side for easier fitting of slings.	Standard Part	
Armrest overhang on inside	Reduces width between armrests by approximately 38 mm (1.5"). Provides extra arm and elbow support if the user is 'pear shaped' or very thin.	Standard Part	
Right hand control (std)	N/A	Standard Part	
Left hand control	N/A	Standard Part	
Carer control	Fitted to backrest out of user reach	Standard Part	
Rise (std)	23° of tilt. The most frequently used rise setting. This provides a gentle, progressive movement.	Standard Part	
Medium rise	14° of tilt. There is not a specific user group for this setting - it is normally used if 'Rise' or 'Flat Rise' is inappropriate.	Standard Part	
Flat rise	2° of tilt and 254 mm (10") of rise, used for users with Muscular Dystrophy who need to lock knees before transferring. (Not available on Small chair)	Standard Part	
Legrest Setting 1.	Standard Position	Standard Part	
Legrest Setting 2.	Standard Position + 25 mm (1")	Standard Part	
Legrest Setting 3.	Standard Position + 50 mm (2")	Standard Part	




Figure 1. Spare parts of Configura® Comfort 1 Chair

WIDTH	711 mm (28")	762 mm (30")	812 mm (32")
LENGTH - UPRIGHT POSITION	940 mm (37")	940 mm (37")	940 mm (37")
LENGTH - RECLINE POSITION	1702 mm (67")	1702 mm (67")	1702 mm (67")
HEIGHT	1092 mm (43")	1143 mm (45")	1143 mm (45")
MASS	65 kg (143 lbs)	68 kg (150 lbs)	71 kg (156 lbs)

SEAT HEIGHT	406 mm (16")	457 mm (18")	457 mm (18")
SEAT HEIGHT WITH (HAKS)	457 mm (18")	508 mm (20")	508 mm (20")
SEAT WIDTH	457 mm (18")	508 mm (20")	558 mm (22")
SEAT DEPTH POSITION	457 mm (18")	508 mm (20")	508 mm (20")
ARMREST HEIGHT	180 mm (7")	180 mm (7")	180 mm (7")
RISE HEIGHT TO FRONT OF CHAIR			
STANDARD	546 mm (21.5")	508 mm (20")	508 mm (20")
MEDIUM	622 mm (24.5")	572 mm (22.5")	572 mm (22.5")
FLAT	NA	711 mm (28")	711 mm (28")
CHAIR TO WALL / FURNITURE	457 mm (18")	483 mm (19")	483 mm (19")
SAFE WORKING LOAD	130 kg	160 kg	160kg
BACK HEIGHT	700 mm	700 mm	700 mm

Clean and care for the fabric on a regular basis and it will maintain its condition and appearance.

The new generation of Microfibre cleaning cloths without soap are extremely effective for general care and maintenance. You can also wipe clean regularly with a damp cloth to remove dust particles and light soiling. It is recommended to vacuum the upholstery once a week.

Most stains can be removed by wiping clean. Remove spillages promptly (within 15 minutes) with cold water. For extreme soiling control use a 5% Sodium Hypochlorite (Bleach) solution.

Do not machine wash or dry clean.

Do not apply other chemicals such as polish, solvents or stain remover.

The integrated Self Supply Overlay comes as standard on all chairs.

This allows for easy integration of existing cushions.

The size of the aperture is 445 mm x 445 mm x 76 mm (17.5" x 17.5" x 3").

This means the height of the seat and armrests remain unchanged when using another cushion.



STEP 1

Remove the backrest cushions.



STEP 2

Unzip the cushion cover to expose the foam.



STEP 3

Pull the foam out carefully, starting at the rear and cut the perforations at the front.



STEP 4

This leaves the cut out for the existing cushions to be inserted. Insert the cushion and re zip the cover.

