

## CONTROL PANEL



Power /  
Mute / Unlock / Reset



Upright



Max. Inflation



Comfort  
Level Setting



Indicator LED



Cycle Time  
(Static Mode & Dynamic Mode)



Cycle Timer LED



Pressure  
Monitoring LED (IPS)

## OPERATING INSTRUCTION



### Power

- Plug the power cord into a wall socket
- Switch on the rocker switch on the side
- Press to power on the unit, and then the green LED will light up
- \* The LED will change to yellow when the power cord is removed without switching off the rocker switch or when there is an abnormal power outage

### Mute

- Press to turn off the audible indicator

### Unlock

- The panel will auto-lock when there is no operation for 2 minutes
- Only and are functional when the panel is locked
- Press and hold for 3 seconds to unlock the panel

### Reset

- Press to re-power and switch back to the default setting (comfort level 0 and dynamic mode in 10 minutes cycle time)



### Comfort Level Setting

- Press to adjust the setting from -2 to 2
- \* The comfort setting scale is only indicative  
If the mattress seems too soft or too rigid, please adjust the comfort setting to conform to each patient's requirements

## Therapy Mode Selection \* The default setting is comfort level 0 and dynamic mode in 10 minutes cycle time



### Cycle Time (Static Mode & Dynamic Mode)


Press to select between Static Mode & Dynamic Mode with different cycle time setting

- Set the cycle time to 0 for Static Mode, in which all air cells maintain constant pressure support and maximize patient's contact area to redistribute pressure
- Set the cycle time to 5, 10, or 15 minutes for Dynamic Mode, in which air cells continuously alternate in an A-B-A-B pattern to relieve pressure and increase blood flow of the patient's tissue
- \* The 3 cells at the head are always excluded from the alternating



## Upright



Upright Mode is used to prevent the patient from bottoming out in an upright position

- Press  to activate Upright Mode
- A green LED will light up to indicate this mode is in operation



## Max. Inflation

The system will rapidly bring the mattress to maximum steady pressure, allowing caregivers to perform nursing procedures

- Press  to activate Max. Inflation
- A yellow LED will light up to indicate the activation of this function
- Press  again to stop this function, or the system will automatically return to the previous setting after 20 minutes



## Cycle Timer LED

The LED light bar shows the countdown of cycle time

\* The light bar will be off in Static Mode and Max. Inflation



## Indicator LED



- An yellow LED will light up if the mattress has abnormal pressure

\* The LED will stay on until the problem has been solved



## Pressure Monitoring LED (IPS)

- **IPS = Intelligent Pressure Sensing** The system monitors the mattress' pressure 24 hours
- The yellow LED indicates this function is in operation

PROBLEM	CONTROL PROCEDURE	POSSIBLE SOLUTION
 The mattress is not inflating even though the power unit is working	<b>A</b> Verify that air flows freely across the tubes and the mattress manifold. Check for any cuts, blockages, or breaks in the air	<b>A</b> If necessary, adjust the position of the tubes or manifold to prevent kinks or twists. Any cuts or rips in the air cells or air hoses, replace them
	<b>B</b> Confirm that the quick coupling is properly connected to the air outlets of the power unit	<b>B</b> Make sure the quick coupling is firmly connected
	<b>C</b> Verify the CPR valve is correctly closed	<b>C</b> Ensure the CPR valve is securely closed
 The patient sinks into the mattress	<b>A</b> Check the comfort level setting on the power unit	<b>A</b> Increase the patient's comfort level setting until the correct support pressure is achieved
	<b>B</b> Check for any abnormal air loss from the mattress	<b>B</b> Replace the abnormally losing air components with an authentic replacement part
	<b>C</b> Check the air filter	<b>C</b> Clean or replace the air filter
	<b>D</b> Verify that the CPR valve is correctly closed	<b>D</b> Firmly close the CPR valve
The power unit cannot power on	<b>A</b> Verify that the power cord set plugs into the proper socket	<b>A</b> Insert the power cord set of power unit into an appropriate socket and turn the power on
	<b>B</b> Verify that the power cord set is securely connected to the power unit	<b>B</b> Insert the power cord set into the power unit and turn the power on
	<b>C</b> Verify that the power cord set is not damaged	<b>C</b> Replace with a functioning power cord set
	<b>D</b> Verify that the fuses are not Blown	<b>D</b> Contact the authorized distributor for technical service
	<b>E</b> The power unit is not responding to the control procedures listed above	<b>E</b> Contact the authorized distributor for technical service