VICAIR® ADJUSTER 02

THE OPTIMAL CUSHION FOR WHEELCHAIR USERS WITH BODY ASYMMETRY

PRESSURE INJURY TREATMENT PLAN
















- Individual solution for asymmetry
- Minimizes sliding and shear forces
- ▲ Optimal adjustment to the body's contours for equal pressure distribution

The Vicair Adjuster O2 consists of 5 individual compartments, designed to contour under loading. Each compartment is designed to influence a specific area of the pelvis, striking a unique balance between preserving skin and underlying tissue at the ischials, coccyx and trochanters, positioning the pelvis and supporting the femurs. For specific needs the filling of each of the 5 compartments can be easily adjusted. Allowing you to change the density of each compartment benefitting users who: (1) need accommodation of a mild to severe fixed pelvic obliquity or correction of a flexible pelvic obliquity; (2) benefit from additional support of the thighs, for instance users with an above knee amputation.

AVAILABLE COVERS

Vicair O2 Top Cover, Comfair Cover, Incotec Cover and Reversible Cover.

PRODUCT PERFORMANCE	
Skin protection / Anti decubitus	
Pressure redistribution	***
Shear force reduction	***
Microclimate control	
Impact & vibration damping	***
Positioning	
Lateral stability	
Frontal stability	***
Sliding reduction	***
Allows manual adjustments	***
Increasing endurance	
Maintenance & Safety	
Durability	***
Reliability	***
Functionality when misused	
Easy to maintain	***
Level of hygiene	***
Functionality & Application	
Comfort	***
Ultra lightweight	***
Allows easy transfers	
Available standard sizes	***
Maximum user weight	

TECHNICAL SPECIFICATIONS								
Thickness	6 and 10 cm							
Average weight	725 grams							
Max. user weight	250 kg							
Machine washable	60°C max.							

Thigh support

CUSHION SIZES												
Standard (Sizes >53 cm thickness 10cm only)												
On request												
Inch		14	15	16	17	18	19	20	21	22	23	widt 24
	см		38		43			50	53	55	58	60
14	35											
15	38											
16	40											
17	43											
18	45											
19	48											
20	50											

21 53

070069-OAD-EN-REV01/2024